



CLVS1 mouse mAb


Catalog No	YP-mAb-07935
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CLVS1 CRALBPL RLBP1L1
Immunogen	Synthesized peptide derived from human CLVS1 AA range: 16-66
Specificity	This antibody detects endogenous levels of CLVS1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.49% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C6orf212L;CRALBPL;RLBP1L1;493342J24Rik;Clv1;CLVS1_HUMAN;CLVS1; Cellular retinaldehyde-binding protein-like;Retinaldehyde-binding protein 1-like 1; clathrin vesicle-associated Sec14 protein 1;CLVS1_MOUSE;clavesin 1; retinaldehyde binding protein 1-like 1;CLVS1 CRALBPL RLBP1L1;Clavesin-1 (Cellular retinaldehyde-binding protein-like) (Retinaldehyde-binding protein 1-like 1)
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Golgi apparatus, trans-Golgi network membrane ; Peripheral membrane protein . Cytoplasmic vesicle, clathrin-coated vesicle . Early endosome membrane ; Peripheral membrane protein .
Tissue Specificity	Expressed mainly in the brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CLVS1 mouse mAb