




F181A mouse mAb

Catalog No	YP-mAb-07938
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM181A C14orf152 HSD-31 HSD31
Immunogen	Synthesized peptide derived from human F181A AA range: 42-92
Specificity	This antibody detects endogenous levels of F181A at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.52% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C14orf152;F181A_HUMAN;FAM181A;Family with sequence similarity 181;HSD-31;HSD31;Protein FAM181A;F181A;FAM181A C14orf152 HSD-31 HSD31
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	
Tissue Specificity	Hippocampus,Ovary,Testis,
Function	similarity:Belongs to the FAM181 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using F181A mouse mAb