




ZN257 mouse mAb

Catalog No	YP-mAb-07999
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF257 BMZF4
Immunogen	Synthesized peptide derived from human ZN257 AA range: 391-441
Specificity	This antibody detects endogenous levels of ZN257 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.113% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZN257;ZNF257 BMZF4;Zinc finger protein 257 (Bone marrow zinc finger 4) (BMZF-4)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Nucleus .
Tissue Specificity	Bone marrow,Testis,
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 12 C2H2-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ZN257 mouse mAb