



CP062 mouse mAb


Catalog No	YP-mAb-08146
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C16orf62 101F10.2
Immunogen	Synthesized peptide derived from human CP062 AA range: 304-354
Specificity	This antibody detects endogenous levels of CP062 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.261% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C16orf62;Chromosome 16 open reading frame 62;CP62_HUMAN; DKFZp313M539;DKFZp434B212;Esophageal cancer associated protein; Esophageal cancer-associated protein;FLJ214;Hypothetical protein LOC572; MGC16824;UPF55 protein C16orf62;CP062;C16orf62 101F10.2;UPF0505 protein C16orf62 (Esophageal cancer-associated protein)
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	Membrane ; Single-pass membrane protein . Endosome . Endosome location is dependent of the association with the CCC and WASH complexes. .
Tissue Specificity	Epithelium,Heart,Pancreas,Testis,Tongue,Uterus,
Function	similarity:Belongs to the UPF0505 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CP062 mouse mAb