




K0753 mouse mAb

Catalog No	YP-mAb-08159
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KIAA0753
Immunogen	Synthesized peptide derived from human K0753 AA range: 778-828
Specificity	This antibody detects endogenous levels of K0753 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.274% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	K753_HUMAN;Kiaa753;Uncharacterized protein KIAA753;K0753;KIAA0753;Uncharacterized protein KIAA0753
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriole . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriolar satellite . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Localization to centrioles and pericentriolar satellites may be mediated by interaction with PCM1. .
Tissue Specificity	Brain,Testis,Uterus,
Function	sequence caution:Translated as Glu.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using K0753 mouse mAb