




EMIL3 mouse mAb

Catalog No	YP-mAb-08198
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EMILIN3 C20orf130 EMILIN5
Immunogen	Synthesized peptide derived from human EMIL3 AA range: 2-52
Specificity	This antibody detects endogenous levels of EMIL3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.313% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EMIL3;EMILIN3 C20orf130 EMILIN5;EMILIN-3 (EMILIN-5) (Elastin microfibril interface-located protein 3) (Elastin microfibril interfacier 3) (Elastin microfibril interface-located protein 5) (Elastin microfibril interfacier 5)
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Secreted, extracellular space, extracellular matrix .
Tissue Specificity	Mesenchymal Stem Cell,PCR rescued clones,Testis,
Function	similarity:Contains 1 EMI domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using EMIL3 mouse mAb