




PDZD4 mouse mAb

Catalog No	YP-mAb-08214
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PDZD4 KIAA1444 PDZK4 PDZRN4L
Immunogen	Synthesized peptide derived from human PDZD4 AA range: 4-54
Specificity	This antibody detects endogenous levels of PDZD4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.329% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA1444;PDZ domain-containing protein 4;PDZ domain-containing RING finger protein 4-like protein;PDZK4;PDZRN4L;PDZD4;PDZD4 KIAA1444 PDZK4 PDZRN4L;PDZ domain-containing protein 4 (PDZ domain-containing RING finger protein 4-like protein)
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Cytoplasm, cell cortex . Mainly localized under the plasma membrane.
Tissue Specificity	Brain-specific. Expressed in fetal and adult brain. Up-regulated in synovial carcinomas.
Function	similarity:Contains 1 PDZ (DHR) domain.,subcellular location:Mainly localized under the plasma membrane.,tissue specificity:Brain-specific. Expressed in fetal and adult brain. Up-regulated in synovial carcinomas.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using PDZD4 mouse mAb