



FTSJD1 mouse mAb

Catalog No	YP-mAb-08219
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FTSJD1 AFT
Immunogen	Synthesized peptide derived from human FTSJ1 AA range: 462-512
Specificity	This antibody detects endogenous levels of FTSJ1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.334% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cap methyltransferase 1;Cap-specific mRNA (nucleoside-2"-O-)-methyltransferase 1;Cap1 2"O-ribose methyltransferase 1;CMTR1;CMTR1_HUMAN;FLJ22156;FtsJ methyltransferase domain containing 2;FtsJ methyltransferase domain-containing protein 2;FTSJD 2;hMTr1;Interferon-stimulated gene 95 kDa protein;ISG95;KIAA82;MTr1;OTTHUMP16324;Uncharacterized protein KIAA82;FTSJD1;FTSJD1 AFT;FTSJ1;FtsJ methyltransferase domain-containing protein 1 (EC 2.1.1.-) (Protein adrift homolog)
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Nucleus . Cytoplasm .
Tissue Specificity	Placenta,Testis,Trachea,
Function	similarity:Belongs to the methyltransferase superfamily. RrmJ family.,
Background	

**matters needing attention**


Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using FTSJD1 mouse mAb