




PSD2 mouse mAb

Catalog No	YP-mAb-08225
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PSD2
Immunogen	Synthesized peptide derived from human PSD2 AA range: 57-107
Specificity	This antibody detects endogenous levels of PSD2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.340% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EFA6C;PSD2;PH and SEC7 domain-containing protein 2 (Pleckstrin homology and SEC7 domain-containing protein 2)
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Cell membrane ; Single-pass membrane protein . Cell projection, ruffle membrane . Cleavage furrow . In interphase associated with the plasma membrane, in particular with membrane ruffling regions. In cells undergoing cytokinesis, transiently found around the ingressing cleavage furrow. Not detected at the midbody ring/Flemming body. .
Tissue Specificity	Amygdala,Brain,
Function	similarity:Belongs to the PSD family.,similarity:Contains 1 PH domain.,similarity:Contains 1 SEC7 domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using PSD2 mouse mAb