




ZBED5 mouse mAb

Catalog No	YP-mAb-08290
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZBED5
Immunogen	Synthesized peptide derived from human ZBED5 AA range: 521-571
Specificity	This antibody detects endogenous levels of ZBED5 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZBED5
Observed Band	
Calculated Molecular Weight	
Cell Pathway	nucleoplasm,cytoplasm,
Tissue Specificity	
Function	miscellaneous:May be derived from an ancient transposon that has lost its ability to translocate.,similarity:Contains 1 BED-type zinc finger.,
Background	This gene is unusual in that its coding sequence is mostly derived from Charlie-like DNA transposon; however, it does not appear to be an active DNA transposon as it is not flanked by terminal inverted repeats. The encoded protein is conserved among the mammalian Laurasiatheria branch. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Jan 2009],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ZBED5 mouse mAb