




EXD1 mouse mAb

Catalog No	YP-mAb-08378
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EXD1 EXDL1
Immunogen	Synthesized peptide derived from human EXD1 AA range: 280-330
Specificity	This antibody detects endogenous levels of EXD1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	piRNA biogenesis protein EXD1;Inactive exonuclease EXD1;Exonuclease 3'-5' domain-like-containing protein 1;Exonuclease 3'-5' domain-containing protein 1; EXDL1;EXD1;EXD1 EXDL1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Component of the meiotic nuage, also named P granule, a germ-cell-specific organelle required to repress transposon activity during meiosis. .
Tissue Specificity	
Function	similarity:Belongs to the EXD1 family.,similarity:Contains 1 3'-5' exonuclease domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using EXD1 mouse mAb