




FSD2 mouse mAb

Catalog No	YP-mAb-08384
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FSD2 SPRYD1
Immunogen	Synthesized peptide derived from human FSD2 AA range: 156-206
Specificity	This antibody detects endogenous levels of FSD2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FSD2;RP11 127F21;SPRYD1;FSD2 SPRYD1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Sarcoplasmic reticulum . Cytoplasm, perinuclear region . In skeletal muscles and striated muscles flanks Z-disks. Partially colocalizes with RYR2 in the sarcoplasmic reticulum. .
Tissue Specificity	
Function	similarity:Contains 1 B30.2/SPRY domain.,similarity:Contains 2 fibronectin type-III domains.,
Background	This gene encodes a protein that belongs to the FN3/SPRY family of proteins. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using FSD2 mouse mAb