




YN3544

ZDH21 mouse mAb

Catalog No	YP-mAb-08387
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZDHHC21
Immunogen	Synthesized peptide derived from human ZDH21 AA range: 177-227
Specificity	This antibody detects endogenous levels of ZDH21 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1 :500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	29kD
Observed Band	
Cell Pathway	Golgi apparatus membrane ; Multi-pass membrane protein . Golgi apparatus, cis-Golgi network membrane ; Multi-pass membrane protein . Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Widely expressed.
Function	catalytic activity:Palmitoyl-CoA + protein-cysteine = S-palmitoyl protein + CoA., domain:The DHHC domain is required for palmitoyltransferase activity., similarity:Belongs to the DHHC palmitoyltransferase family.,similarity:Contains 1 DHHC-type zinc finger.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ZDH21 mouse mAb