




RASL3 mouse mAb

Catalog No	YP-mAb-08398
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RASAL3
Immunogen	Synthesized peptide derived from human RASL3 AA range: 297-347
Specificity	This antibody detects endogenous levels of RASL3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RASL3;RASAL3
Observed Band	111kD
Calculated Molecular Weight	111kD
Cell Pathway	Cytoplasm . Cytoplasm, cell cortex .
Tissue Specificity	Predominantly expressed in cells of hematopoietic lineages.
Function	function:May functions as a Ras GTPase-activating protein.,sequence caution:Translated as Arg.,sequence caution:Translation N-terminally shortened., similarity:Contains 1 C2 domain.,similarity:Contains 1 PH domain., similarity:Contains 1 Ras-GAP domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using RASL3 mouse mAb