



GABP2 mouse mAb


Catalog No	YP-mAb-08434
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GABPB2
Immunogen	Synthesized peptide derived from human GABP2 AA range: 314-364
Specificity	This antibody detects endogenous levels of GABP2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BABPB2;E4TF1;E4TF1-47;E4TF1-53;E4TF1B;GA binding protein transcription factor beta subunit 1 transcript variant gamma 2;GA binding protein transcription factor;beta subunit 1;beta subunit 2;GA-binding protein subunit beta-1;GABP subunit beta-1;GABP subunit beta-2;GABP1_HUMAN;GABPB;GABPB-1;GABPB-2;GABPB1;GABPB2;NRF2B1;NRF2B2;Nuclear respiratory factor 2; Transcription factor E4TF1-47;Transcription factor E4TF1-53;GABP2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:May function as transcription factor caMABLE of interacting with purine rich repeats (GA repeats).,similarity:Contains 5 ANK repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using GABP2 mouse mAb