




CC020 mouse mAb

Catalog No	YP-mAb-08442
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C3orf20
Immunogen	Synthesized peptide derived from human CC020 AA range: 272-322
Specificity	This antibody detects endogenous levels of CC020 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C3orf2;CC2_HUMAN;Chromosome 3 open reading frame 2;DKFZP434N1817; Uncharacterized protein C3orf2;CC020;C3orf20
Observed Band	99kD
Calculated Molecular Weight	99kD
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CC020 mouse mAb