



K1614 mouse mAb


Catalog No	YP-mAb-08461
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KIAA1614
Immunogen	Synthesized peptide derived from human K1614 AA range: 775-825
Specificity	This antibody detects endogenous levels of K1614 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	63441H18Rik;Coiled coil domain containing protein MTMR15;DKFZp451H236; DKFZp686K16147;FAN1;FAN1_HUMAN;FANCD2 and FANCI associated nuclease 1;FANCD2/FANCI associated nuclease 1;FANCD2/FANCI-associated nuclease 1;Fanconi associated nuclease 1;Fanconi-associated nuclease 1; hFAN1;KIAA118;KMIN;MTMR 15;MTMR15;Myotubularin related protein 15; Myotubularin-related protein 15;K1614;KIAA1614
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using K1614 mouse mAb