




T106B mouse mAb

Catalog No	YP-mAb-08593
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMEM106B
Immunogen	Synthesized peptide derived from human T106B AA range: 72-122
Specificity	This antibody detects endogenous levels of T106B at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	23136D22Rik;583455K21Rik;643519M21Rik;AI428776;AI661344;FLJ11273;LRRGT11;MGC33727;MGC94135;T16B_HUMAN;Tmem16b;Transmembrane protein 16B;T106B;TMEM106B
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Late endosome membrane; Single-pass type II membrane protein. Lysosome membrane; Single-pass type II membrane protein. Membrane ; Lipid-anchor .
Tissue Specificity	Expressed in frontal cortex.
Function	similarity:Belongs to the TMEM106 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using T106B mouse mAb