




MA2B2 mouse mAb

Catalog No	YP-mAb-08652
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MAN2B2 KIAA0935
Immunogen	Synthesized peptide derived from human MA2B2 AA range: 906-956
Specificity	This antibody detects endogenous levels of MA2B2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MA2B2;MAN2B2 KIAA0935
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	catalytic activity:Hydrolysis of terminal, non-reducing alpha-D-mannose residues in alpha-D-mannosides.,cofactor:Bounds 1 zinc ion per subunit.,similarity:Belongs to the glycosyl hydrolase 38 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using MA2B2 mouse mAb