



BAHD1 mouse mAb


Catalog No	YP-mAb-08718
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BAHD1 KIAA0945
Immunogen	Synthesized peptide derived from human BAHD1 AA range: 214-264
Specificity	This antibody detects endogenous levels of BAHD1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1 :500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AL22997;AW541238;BAH domain containing protein 1;BAH domain containing protein1;BAHDH 1;Bromo adjacent homology domain containing 1;Bromo adjacent homology domain containing1;Bromo adjacent homology domain-containing 1 protein;Gm117;KIAA945;MGC117747;mKIAA945;RGD1563127;BAHD1;BAHD1 KIAA0945
Observed Band	86kD
Calculated Molecular Weight	86kD
Cell Pathway	Nucleus . Chromosome . Localizes to heterochromatin and inactive X chromosome. Colocalizes with histone H3 trimethylated at 'Lys-27' (H3K27me3).
Tissue Specificity	
Function	similarity:Contains 1 BAH domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from VEC cells,
primary antibody was diluted at 1:1000, 4° over night