




ARP5L mouse mAb

Catalog No	YP-mAb-08803
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ARPC5L
Immunogen	Synthesized peptide derived from human ARP5L AA range: 4-54
Specificity	This antibody detects endogenous levels of ARP5L at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ARP5L;ARPC5L
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton .
Tissue Specificity	
Function	function:May function as component of the Arp2/3 complex which is involved in regulation of actin polymerization and together with an activating nucleation-promoting factor (NPF) mediates the formation of branched actin networks., similarity:Belongs to the ARPC5 family.,subunit:May be a component of the Arp2/3 complex in which it may replace ARPC5.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ARP5L mouse mAb