




KANK4 mouse mAb

Catalog No	YP-mAb-08915
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KANK4 ANKRD38
Immunogen	Synthesized peptide derived from human KANK4 AA range: 448-498
Specificity	This antibody detects endogenous levels of KANK4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANKRD38;Ankyrin repeat domain 38;Ankyrin repeat domain-containing protein 38;dJ178M7.1;Kank4;KANK4_HUMAN;Kidney ankyrin repeat containing protein 4;KN motif and ankyrin repeat domain-containing protein 4;KN motif and ankyrin repeat domains 4;RP5-1155K23.5;KANK4 ANKRD38
Observed Band	
Calculated Molecular Weight	109kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Strongly expressed in colon, liver, lung, skeletal muscle and kidney.
Function	similarity:Contains 5 ANK repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using KANK4 mouse mAb