




HSDL1 mouse mAb

Catalog No	YP-mAb-08933
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HSDL1
Immunogen	Synthesized peptide derived from human HSDL1 AA range: 83-133
Specificity	This antibody detects endogenous levels of HSDL1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HSDL1;SDR12C3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion .
Tissue Specificity	Highly expressed in testis and ovary. Also detected in thyroid, spinal cord, adrenal gland, heart, placenta, skeletal muscle, small intestine, colon, spleen, prostate and pancreas.
Function	caution:Although related to the SDR family, lacks the conserved active Tyr residue in position 218 which is replaced by a Phe, suggesting that it may lack oxidoreductase activity.,similarity:Belongs to the short-chain dehydrogenases/reductases (SDR) family. 17-beta-HSD 3 subfamily.,tissue specificity:Highly expressed in testis and ovary. Also detected in thyroid, spinal cord, adrenal gland, heart, placenta, skeletal muscle and pancreas.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using HSDL1 mouse mAb