




DIP2A mouse mAb

Catalog No	YP-mAb-08994
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DIP2A C21orf106 DIP2 KIAA0184
Immunogen	Synthesized peptide derived from human DIP2A AA range: 1080-1130
Specificity	This antibody detects endogenous levels of DIP2A at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DIP2A;DIP2A C21orf106 DIP2 KIAA0184
Observed Band	173kD
Calculated Molecular Weight	173kD
Cell Pathway	Cell membrane ; Peripheral membrane protein . Mitochondrion . Cell projection, dendritic spine .
Tissue Specificity	Low expression in all tissues tested.
Function	function:May provide positional cues for axon pathfinding and patterning in the central nervous system.,similarity:Belongs to the DIP2 family.,tissue specificity:Low expression in all tissues tested.,
Background	The protein encoded by this gene may be involved in axon patterning in the central nervous system. This gene is not highly expressed. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2009],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using DIP2A mouse mAb