




LACC1 mouse mAb

Catalog No	YP-mAb-09053
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LACC1 C13orf31
Immunogen	Synthesized peptide derived from human LACC1 AA range: 4-54
Specificity	This antibody detects endogenous levels of LACC1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	chromosome 13 open reading frame 31;LACC1;LACC1 C13orf31
Observed Band	
Calculated Molecular Weight	47kD
Cell Pathway	Cytoplasm . Nucleus . Endoplasmic reticulum . Peroxisome . Upon stimulation of the pattern-recognition receptor (PRR) NOD2, localizes to the endoplasmic reticulum. .
Tissue Specificity	Ubiquitously expressed, with higher expression levels in immune-related tissues such as lymph nodes and spleen (PubMed:27959965). Expressed in both intestinal and peripheral myeloid-derived cells (PubMed:28593945).
Function	similarity:Belongs to the UPF0124 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using LACC1 mouse mAb