



SVIP mouse mAb

Catalog No	YP-mAb-09181
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SVIP
Immunogen	Synthesized peptide derived from human SVIP AA range: 26-76
Specificity	This antibody detects endogenous levels of human SVIP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Small VCP/p97-interacting protein;SVIP
Observed Band	12 kDa
Calculated Molecular Weight	8kDa
Cell Pathway	Smooth endoplasmic reticulum membrane ; Peripheral membrane protein . Golgi apparatus membrane ; Peripheral membrane protein . Cell membrane ; Peripheral membrane protein . Membrane ; Lipid-anchor .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

