



KYNU mouse mAb

Catalog No	YP-mAb-09211
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KYNU
Immunogen	Synthesized peptide derived from human KYNU AA range: 120-170
Specificity	This antibody detects endogenous levels of human KYNU
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	465aa,52 kDa; 307aa,35 kDa
Observed Band	45-52 kDa
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	Expressed in all tissues tested (heart, brain placenta, lung, liver, skeletal muscle, kidney and pancreas). Highest levels found in placenta, liver and lung. Expressed in all brain regions.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

