



TPGS1 mouse mAb

Catalog No	YP-mAb-09644
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TPGS1 C19orf20
Immunogen	Synthesized peptide derived from human TPGS1 AA range: 158-208
Specificity	This antibody detects endogenous levels of human TPGS1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	32kD
Observed Band	29 kDa
Cell Pathway	Cytoplasm, cytoskeleton, cilium axoneme . Cytoplasm, cytoskeleton, flagellum axoneme . Cytoplasm, cytoskeleton, cilium basal body . Cytoplasm, cytoskeleton, flagellum basal body . Cell projection, axon . Cell projection, dendrite . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriolar satellite . Associated with microtubules from neurites, centrosomes, basal bodies and axonemes. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images