



SULF2 mouse mAb

Catalog No	YP-mAb-09758
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SULF2 KIAA1247 UNQ559/PRO1120
Immunogen	Synthesized peptide derived from human SULF2 AA range: 690-740
Specificity	This antibody detects endogenous levels of human SULF2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Arylsulfatase;EC:3.1.6.1;EC:3.1.6.14;Extracellular sulfatase Sulf 2;Extracellular sulfatase Sulf-2;SULF2;SULF2 KIAA1247 UNQ559/PRO1120;Extracellular sulfatase Sulf-2 (hSulf-2) (EC 3.1.6.-)
Observed Band	50 kDa, 130-150 kDa
Calculated Molecular Weight	870 aa, 96 kDa
Cell Pathway	Endoplasmic reticulum . Golgi apparatus, Golgi stack . Cell surface . Also localized on the cell surface. .
Tissue Specificity	Expressed at highest levels in the ovary, skeletal muscle, stomach, brain, uterus, heart, kidney and placenta.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

