



D19L2 mouse mAb

Catalog No	YP-mAb-09913
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DPY19L2 UNQ3127/PRO10284
Immunogen	Synthesized peptide derived from human D19L2 AA range: 246-296
Specificity	This antibody detects endogenous levels of human D19L2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	D19L2;DPY19L2 UNQ3127/PRO10284;Protein dpy-19 homolog 2 (Dpy-19-like protein 2)
Observed Band	83kD
Calculated Molecular Weight	83kD
Cell Pathway	Nucleus inner membrane ; Multi-pass membrane protein . Colocalizes with DPY19L2 at the inner nuclear membrane. .
Tissue Specificity	Widely expressed with high expression in testis. Not detectable in ejaculated sperm (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

