



RN212 mouse mAb

Catalog No	YP-mAb-09932
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RNF212
Immunogen	Synthesized peptide derived from human RN212 AA range: 181-231
Specificity	This antibody detects endogenous levels of human RN212
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	33kD
Synonyms	RN212;RNF212;RING finger protein 212
Observed Band	
Cell Pathway	Nucleus . Chromosome . Associates to the synaptonemal complex. Localizes to a minority of double-strand breaks (DSBs) sites. Marks crossover sites during midpachynema (By similarity). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

