



CERS6 mouse mAb

Catalog No	YP-mAb-09946
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CERS6 LASS6
Immunogen	Synthesized peptide derived from human CERS6 AA range: 73-123
Specificity	This antibody detects endogenous levels of human CERS6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CERS6;CERS6 LASS6;Ceramide synthase 6 (CerS5) (LAG1 longevity assurance homolog 6)
Observed Band	
Calculated Molecular Weight	42kD
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

