



SYLM mouse mAb

Catalog No	YP-mAb-09984
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LARS2 KIAA0028
Immunogen	Synthesized peptide derived from human SYLM AA range: 832-882
Specificity	This antibody detects endogenous levels of human SYLM
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	102 kDa
Observed Band	102kD
Cell Pathway	Mitochondrion matrix .
Tissue Specificity	Ubiquitously expressed, but highest expression in tissues with high metabolic rates, such as skeletal muscle, heart, and kidney.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

