



ANKH1 mouse mAb

Catalog No	YP-mAb-10046
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ANKHD1 KIAA1085 MASK VBARP PP2500
Immunogen	Synthesized peptide derived from human ANKH1 AA range: 435-485
Specificity	This antibody detects endogenous levels of human ANKH1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANKH1;ANKHD1 KIAA1085 MASK VBARP PP2500;Ankyrin repeat and KH domain-containing protein 1 (HIV-1 Vpr-binding ankyrin repeat protein) (Multiple ankyrin repeats single KH domain) (hMASK)
Observed Band	
Calculated Molecular Weight	280kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Ubiquitous with high expression in cervix, spleen and brain. Expressed in hematopoietic cells with increased expression in leukemia cells. Isoform 2 is highly expressed in spleen with almost no expression in muscle and brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

