



FOXK2 mouse mAb

Catalog No	YP-mAb-10102
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FOXK2 ILF ILF1
Immunogen	Synthesized peptide derived from human FOXK2 AA range: 331-381
Specificity	This antibody detects endogenous levels of human FOXK2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	forkhead box K2;Forkhead box protein K2;ILF;ILF 1;ILF1;FOXK2;FOXK2 ILF ILF1;Forkhead box protein K2 (Cellular transcription factor ILF-1) (FOXK1) (Interleukin enhancer-binding factor 1)
Observed Band	69-79 kDa
Calculated Molecular Weight	69 kDa
Cell Pathway	Nucleus . Cytoplasm .
Tissue Specificity	Expressed in both lymphoid and non-lymphoid cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

