



MPPD2 mouse mAb

Catalog No	YP-mAb-10179
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MPPED2 C11orf8 FAM1B
Immunogen	Synthesized peptide derived from human MPPD2 AA range: 143-193
Specificity	This antibody detects endogenous levels of human MPPD2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MPPD2;MPPED2 C11orf8 FAM1B;Metallophosphoesterase MPPED2 (EC 3.1.-.-) (Fetal brain protein 239) (239FB) (Metallophosphoesterase domain-containing protein 2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed predominantly in fetal brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using MPPD2 mouse mAb

