



ZIC5 mouse mAb

Catalog No	YP-mAb-10196
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZIC5
Immunogen	Synthesized peptide derived from human ZIC5 AA range: 566-616
Specificity	This antibody detects endogenous levels of human ZIC5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Odd paired homolog Drosophila;Zic 1;ZIC;Zic family member 1 (odd-paired Drosophila homolog);Zic family member 1;Zic protein member 1;zic1; ZIC1_HUMAN;Zinc finger protein 201;Zinc finger protein of the cerebellum 1; Zinc finger protein ZIC 1;Zinc finger protein ZIC1;ZNF 201;ZNF201;HPE 5;HPE5; Zic 2;Zic family member 2 (odd paired Drosophila homolog);Zic family member 2; ZIC2;ZIC2_HUMAN;Zinc finger protein of the cerebellum 2;Zinc finger protein ZIC 2;Zinc finger protein Zic2;Heterotaxy 1;HTX;HTX1;ZIC 3;Zic family member 3 (odd paired Drosophila homolog heterotaxy 1);Zic family member 3;Zic3; ZIC3_HUMAN;Zinc finger protein 203;Zinc finger protein of the cerebellum 3; Zinc finger protein ZIC 3 (Zinc finger protein of the cerebellum 3);Zinc finger protein ZIC 3;ZNF203;ZIC5;Zinc finger protein ZIC 5 (Zinc finger protein of the cerebellum 5)
Observed Band	51kDa
Calculated Molecular Weight	48kD,55kD,51kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

