



MFAP2 mouse mAb

Catalog No	YP-mAb-10224
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MFAP2 MAGP1
Immunogen	Synthesized peptide derived from human MFAP2 AA range: 101-151
Specificity	This antibody detects endogenous levels of human MFAP2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MFAP2;MFAP2 MAGP1;Microfibrillar-associated protein 2 (MFAP-2) (Microfibril-associated glycoprotein 1) (MAGP) (MAGP-1)
Observed Band	20kD
Calculated Molecular Weight	20kD
Cell Pathway	Secreted, extracellular space, extracellular matrix.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

