




Tubulin alpha mouse mAb(Caenorhabditis elegans)

Catalog No	YP-mAb-10299
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Purified recombinant Full length Tubulin alpha (Caenorhabditis elegans)
Specificity	This antibody detects endogenous levels of Tubulin alpha (Caenorhabditis elegans)
Formulation	Liquid in PBS containing 50% glycerol and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TUBA1B;Tubulin alpha 1B chain;Alpha-tubulin ubiquitous;EC:3.6.5.-;Tubulin alpha;Alpha Tubulin
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using Tubulin alpha mouse mAb(Caenorhabditis elegans)