



IL-31 Monoclonal Antibody

Catalog No	YP-mAb-10637
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	IL31
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human IL31. AA range:21-70
Specificity	The antibody detects endogenous IL-31
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IL 31;IL31;IL-31;interleukin 31;Interleukin-31 (IL-31)
Observed Band	65kD
Calculated Molecular Weight	65kD
Cell Pathway	Secreted .
Tissue Specificity	Detected at low levels in testis, bone marrow, skeletal muscle, kidney, colon, thymus, small intestine and trachea.
Function	function:Activates STAT3 and possibly STAT1 and STAT5 through the IL31 heterodimeric receptor composed of IL31RA and OSMR. IL31 may function in skin immunity.,induction:Up-regulated in skin homing T-cells of patients with atopic dermatitis (AD) and in activated circulating T-cells. Up-regulated in lesional biopsies of patients with allergic contact dermatitis (ACD).,tissue specificity:Detected at low levels in testis, bone marrow, skeletal muscle, kidney, colon, thymus, small intestine and trachea.,
Background	IL31, which is made principally by activated Th2-type T cells, interacts with a heterodimeric receptor consisting of IL31RA (MIM 609510) and OSMR (MIM 601743) that is constitutively expressed on epithelial cells and keratinocytes.



IL31 may be involved in the promotion of allergic skin disorders and in regulating other allergic diseases, such as asthma (Dillon et al., 2004 [PubMed 15184896]) .[supplied by OMIM, Mar 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using IL-31 Monoclonal Antibody