



CTXL Monoclonal Antibody


Catalog No	YP-mAb-10755
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VSIG2 CTH CTXL UNQ2770/PRO7154
Immunogen	Synthesized peptide derived from human CTXL. at AA range: 260-290
Specificity	This antibody detects endogenous levels of CTXL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	221413P1Rik;CT like protein;CTH;CTXL;VSIG2;VSIG2 CTH CTXL UNQ2770/PRO7154;V-set and immunoglobulin domain-containing protein 2 (Cortical thymocyte-like protein) (CT-like protein)
Observed Band	34kD
Calculated Molecular Weight	34kD
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Highly expressed in stomach, colon, prostate, trachea and thyroid glands and weakly in bladder and lung.
Function	similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain., similarity:Contains 1 Ig-like V-type (immunoglobulin-like) domain.,tissue specificity:Highly expressed in stomach, colon, prostate, trachea and thyroid glands and weakly in bladder and lung.,
Background	similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain., similarity:Contains 1 Ig-like V-type (immunoglobulin-like) domain.,tissue specificity:Highly expressed in stomach, colon, prostate, trachea and thyroid glands and weakly in bladder and lung.,

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using CTXL Monoclonal Antibody