



SH3BGRL3 Monoclonal Antibody


Catalog No	YP-mAb-10910
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SH3BGRL3 P1725
Immunogen	Synthesized peptide derived from human SH3BGRL3 AA range: 34-84
Specificity	This antibody detects endogenous levels of human SH3BGRL3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SH3 domain binding protein 1;SH3BGRL3;SH3BP 1;TIP B1;SH3L3;SH3BGRL3 P1725;SH3 domain-binding glutamic acid-rich-like protein 3 (SH3 domain-binding protein 1;SH3BP-1)
Observed Band	
Calculated Molecular Weight	10kD
Cell Pathway	Cytoplasm. Nucleus.
Tissue Specificity	Ubiquitous.
Function	function:Could act as a modulator of glutaredoxin biological activity., similarity:Belongs to the SH3BGR family.,similarity:Contains 1 glutaredoxin domain.,tissue specificity:Ubiquitous.,
Background	function:Could act as a modulator of glutaredoxin biological activity., similarity:Belongs to the SH3BGR family.,similarity:Contains 1 glutaredoxin domain.,tissue specificity:Ubiquitous.,
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using SH3BGRL3 Monoclonal Antibody