




LIGO4 mouse mAb

Catalog No	YP-mAb-10962
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LINGO4 LRRN6D UNQ9248/PRO34002
Immunogen	Synthesized peptide derived from human LIGO4 AA range: 374-424
Specificity	This antibody detects endogenous levels of LIGO4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LIGO4;LINGO4 LRRN6D UNQ9248/PRO34002
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain., similarity:Contains 13 LRR (leucine-rich) repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using LIGO4 mouse mAb