




LCE1D mouse mAb

Catalog No	YP-mAb-11005
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LCE1D LEP4
Immunogen	Synthesized peptide derived from human LCE1D AA range: 56-106
Specificity	This antibody detects endogenous levels of LCE1D at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LCE1D;LCE1D LEP4
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cornified envelope,cytoplasm,perinuclear region of cytoplasm,
Tissue Specificity	Skin-specific. Expression was readily detected in adult trunk skin, adult arm skin, fetal skin, penal skin, vulva, esophagus and tongue. Not expressed in the cervix, rectum, lung, colon, or placenta.
Function	function:Precursors of the cornified envelope of the stratum corneum., induction:By UVB.,miscellaneous:Belongs to the LCE cluster present on 1q21., similarity:Belongs to the LCE family.,tissue specificity:Skin-specific. Expression was readily detected in adult trunk skin, adult arm skin, fetal skin, penal skin, vulva, esophagus and tongue. Not expressed in the cervix, rectum, lung, colon, or placenta.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using LCE1D mouse mAb