




FA13A mouse mAb

Catalog No	YP-mAb-11007
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM13A FAM13A1 KIAA0914
Immunogen	Synthesized peptide derived from human FA13A AA range: 38-88
Specificity	This antibody detects endogenous levels of FA13A at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA13A;FAM13A FAM13A1 KIAA0914
Observed Band	
Calculated Molecular Weight	113kD
Cell Pathway	cytosol,
Tissue Specificity	Isoform 1 is widely expressed, with highest expression in skeletal muscle, thymus, brain and lung. Isoform 3 is less abundant than isoform 1 and predominantly expressed in kidney, pancreas, liver, lung and thymus.
Function	miscellaneous:The sequence shown here is derived from an EMBL/GenBank/ DDBJ third party annotation (TPA) entry.,similarity:Belongs to the FAM13 family., similarity:Contains 1 Rho-GAP domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using FA13A mouse mAb