




SIGL1 mouse mAb

Catalog No	YP-mAb-11030
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SIGLECL1 C19orf75
Immunogen	Synthesized peptide derived from human SIGL1 AA range: 82-132
Specificity	This antibody detects endogenous levels of SIGL1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SIGL1;SIGLECL1 C19orf75
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using SIGL1 mouse mAb