




FA86A mouse mAb

Catalog No	YP-mAb-11048
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM86A SB153
Immunogen	Synthesized peptide derived from human FA86A AA range: 268-318
Specificity	This antibody detects endogenous levels of FA86A at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA86A;FAM86A SB153
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	similarity:Belongs to the FAM86 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using FA86A mouse mAb