




KCD15 mouse mAb

Catalog No	YP-mAb-11051
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KCTD15
Immunogen	Synthesized peptide derived from human KCD15 AA range: 201-251
Specificity	This antibody detects endogenous levels of KCD15 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KCD15;KCTD15
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . In the brain, localizes to the arcuate hypothalamic nucleus, the ventromedial hypothalamic nucleus and the accumbens nucleus of the ventral striatum. .
Tissue Specificity	
Function	similarity:Contains 1 BTB (POZ) domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using KCD15 mouse mAb