



TM165 mouse mAb


Catalog No	YP-mAb-11059
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMEM165 TPARL
Immunogen	Synthesized peptide derived from human TM165 AA range: 68-118
Specificity	This antibody detects endogenous levels of TM165 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM165;TMEM165 TPARL
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus membrane ; Multi-pass membrane protein . Golgi apparatus, trans-Golgi network membrane. Lysosome membrane. Early endosome membrane. Late endosome membrane.
Tissue Specificity	Ubiquitously expressed.
Function	similarity:Belongs to the UPF0016 family.,
Background	This gene encodes a predicted transmembrane protein with a perinuclear Golgi-like distribution in fibroblasts. Mutations in this gene are associated with the autosomal recessive disorder congenital disorder of glycosylation, type IIk. Knockdown of this gene's expression causes decreased sialylation in HEK cells and suggests this gene plays a role in terminal Golgi glycosylation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2012],

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using TM165 mouse mAb