




CH082 mouse mAb

Catalog No	YP-mAb-11080
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C8orf82
Immunogen	Synthesized peptide derived from human CH082 AA range: 163-213
Specificity	This antibody detects endogenous levels of CH082 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CH082;C8orf82
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity:Belongs to the UPF0598 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CH082 mouse mAb